

pGEM®-5Zf(–) Vector



Technical Bulletin No. 068

INSTRUCTIONS FOR USE OF PRODUCT P2351. PLEASE DISCARD PREVIOUS VERSIONS.

All technical literature is available on the Internet at www.promega.com

Please visit the web site to verify that you are using the most current version of this Technical Bulletin.

I. Description	1
II. Product Components	1
III. pGEM®-5Zf(–) Vector Multiple Cloning Region and Circle Map	2
IV. pGEM®-5Zf(–) Vector Restriction Sites	4
V. pGEM®-5Zf(–) Vector Sequence	5
VI. Related Products	7
VII. Reference	8

I. Description

The pGEM®-5Zf(–) Vector^(a) is a derivative of the pGEM®-5Zf(+) Vector^(a) and differs from it by the orientation of the origin of replication of filamentous phage f1. The plasmid serves as a standard cloning vector, as a template for in vitro transcription, and can be used for the production of circular ssDNA. The plasmid contains T7 and SP6 RNA polymerase promoters flanking a multiple cloning region within the α -peptide coding region of β -galactosidase (1). Insertional inactivation of the α -peptide allows recombinant clones to be directly identified by color screening on indicator plates. The multiple cloning region is unique and includes restriction sites for *Apa* I, *Aat* II, *Sph* I, *Nco* I, *Sac* II, *Eco*R V, *Spe* I, *Not* I, *Pst* I, *Sal* I, *Nde* I, *Sac* I, *Bst*X I and *Nsi* I. This arrangement is designed specifically for generating unidirectional deletions with Promega's Erase-a-Base® System. The polylinker contains restriction enzyme sites that produce 5' overhangs or blunt ends (sensitive to Exonuclease III) flanked on both sides by blocks of restriction sites that generate 3' overhangs (resistant to Exonuclease III).

For induction of ssDNA, bacterial cells containing pGEM®-5Zf(–) recombinants are infected with an appropriate helper phage. The plasmid then enters the f1 replication mode, and the resulting ssDNA is exported from the cell as an encapsidated virus-like particle. The sequence of the ssDNA rescued upon infection with helper phage is identical to the sequence shown in Figure 1. The exported ssDNA can be used for mutagenesis in vitro or can be sequenced using Promega's SP6 Promoter Primer (Cat.# Q5011) and pUC/M13 Reverse Primer (Cat.# Q5401, Q5421).

The sequences of Promega vectors are available online at www.promega.com/vectors/ and are also available from the GenBank® database.

II. Product Components

Product	Size	Cat.#
pGEM®-5Zf(–) Vector	20 μ g	P2351

The pGEM®-5Zf(–) Vector is provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain the vector and are not competent.

Storage Conditions: Store the pGEM®-5Zf(–) Vector at –20°C and the glycerol stock of JM109 cells at –70°C.



AF9TB068 0500TB068

III. pGEM®-5Zf(-) Vector Multiple Cloning Region and Circle Map




Figure 1. pGEM®-5Zf(-) Vector promoter and multiple cloning site sequence. The sequence shown corresponds to RNA synthesized by T7 RNA polymerase and is complementary to RNA synthesized by SP6 RNA polymerase. The strand shown is the same as the ssDNA strand produced by this vector.




Figure 2. pGEM®-5Zf(-) Vector circle map and sequence reference points.

pGEM®-5Zf(-) Vector sequence reference points.

T7 RNA polymerase transcription initiation site	1
SP6 RNA polymerase transcription initiation site	126
T7 RNA polymerase promoter (-17 to +3)	2985-3
SP6 RNA polymerase promoter (-17 to +3)	124-143
multiple cloning region	10-113
<i>lacZ</i> start codon	165
<i>lac</i> operon sequences	2822-2982; 151-380
<i>lac</i> operator	185-201
β -lactamase (<i>Amp^r</i>) coding region	1322-2182
phage f1 region	2366-2821
binding site of pUC/M13 Forward Sequencing Primer	2942-2958
binding site of pUC/M13 Reverse Sequencing Primer	161-177

Specialized applications of the pGEM®-5Zf(-) Vector.

- Used with the Erase-a-Base® System.
- ssDNA production.
- Blue/white screening for recombinants.
- Transcription in vitro from dual opposed promoters (For protocol information, please request Promega's *Riboprobe® in vitro Transcription Systems*^(b,c) *Technical Manual*, #TM016.)
- Translation in vitro (For protocol information, please request Promega's *TnT® Quick Coupled Transcription/Translation System*^(b,c,d,e,f) *Technical Manual*, #TM045.)



Promega's
pGEM®-5Zf(-) and
pGEM®-5Zf(+) Vectors
are identical except for
the orientation of the f1
origin.



Use the
SP6 or pUC/M13 Reverse
Primer to sequence
ssDNA produced by the
pGEM®-5Zf(-) Vector.

Note: All Promega
technical literature is
available on the Internet
at www.promega.com.



IV. pGEM®-5Zf(–) Vector Restriction Sites

The following restriction enzyme tables were constructed using DNASTAR® sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. Vector sequences are also available in the GenBank® database (GenBank®/EMBL Accession Number X65309) and on the Internet at www.promega.com/vectors/.

Table 1. Restriction Enzymes That Cut the pGEM®-5Zf(–) Vector Between 1 and 5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Aaf II	1	20	<i>Fsp</i> I	2	1617, 2841
Acc I	1	76	Hae II	4	380, 750, 2441, 2449
Acy I	2	17, 1932	<i>Hga</i> I	4	613, 1191, 1921, 2374
<i>Afl</i> III	2	99, 502	Hinc II	1	77
Alw 26 I	2	1456, 2232	<i>Hind</i> II	1	77
Alw 44 I	2	816, 2062	Hsp 92 I	2	17, 1932
<i>Alw</i> N I	1	918	<i>Mae</i> I	5	56, 997, 1250, 1585, 2443
Apa I	1	14	Mlu I	1	99
<i>Asp</i> H I	4	94, 820, 1981, 2066	Nae I	1	2493
Ava II	2	1533, 1755	Nci I	4	30, 882, 1578, 1929
Ban I	3	246, 1343, 2555	Nco I	1	37
Ban II	3	14, 94, 2525	Nde I	1	82
Bbu I	1	26	Ngo M IV	1	2491
Bgl I	3	39, 1515, 2834	Not I	1	62
<i>Bsa</i> I	1	1456	Nsi I	1	112
<i>Bsa</i> A I	1	2596	<i>Nsp</i> I	2	26, 506
<i>Bsa</i> H I	2	17, 1932	<i>Pae</i> R7 I	1	2560
<i>Bsa</i> J I	5	37, 43, 241, 662, 2937	<i>Ppu</i> 10 I	1	108
<i>Bsp</i> 120 I	1	10	Pst I	1	73
<i>Bsp</i> H I	2	1222, 2230	Pvu I	2	1765, 2862
<i>Bsp</i> M I	1	62	Pvu II	2	326, 2891
<i>Bss</i> S I	2	675, 2059	Rsa I	1	1875
Bst O I	5	242, 530, 651, 664, 2938	Sac I	1	94
Bst X I	1	103	Sac II	1	46
Bst Z I	2	31, 62	Sal I	1	75
<i>Cfr</i> 10 I	2	1475, 2491	Sca I	1	1875
Dde I	4	777, 1186, 1352, 1892	Sfi I	1	39
Dra I	3	1261, 1280, 1972	Sin I	2	1533, 1755
<i>Dra</i> III	1	2599	Spe I	1	55
<i>Drd</i> I	2	610, 2643	Sph I	1	26
<i>Dsa</i> I	2	37, 43	<i>Sse</i> 8387 I	1	73
<i>Eag</i> I	2	31, 62	Ssp I	2	2199, 2804
<i>Ear</i> I	3	386, 2190, 2879	Sty I	1	37
Ec /HK I	1	1395	Taq I	4	76, 602, 2046, 2561
Eco 52 I	2	31, 62	<i>Tfi</i> I	2	337, 477
Eco CR I	1	92	Vsp I	3	273, 332, 1567
Eco R V	1	51	Xmn I	1	1994
Fok I	5	119, 1361, 1542, 1829, 2917			

Note: The enzymes listed in boldface type are available from Promega.

Table 2. Restriction Enzymes That Do Not Cut the pGEM®-5Zf(–) Vector.

Acc III	Bcl I	Bsu36 I	Hpa I	<i>PshA I</i>	Xba I
AccB7 I	Bgl II	Cla I	I-Ppo I	<i>Psp5 II</i>	<i>Xcm I</i>
Acc65 I	<i>Blp I</i>	Csp I	<i>Kas I</i>	<i>PspA I</i>	Xho I
<i>Afl II</i>	<i>Bpu1102 I</i>	Csp45 I	Kpn I	<i>Rsr II</i>	Xma I
Age I	<i>BsaB I</i>	<i>Dra II</i>	Nar I	Sgf I(g)	
<i>Asc I</i>	BsaM I	Eco47 III	Nhe I	<i>SgrA I</i>	
Ava I	<i>Bsm I</i>	<i>Eco72 I</i>	Nru I	Sma I	
<i>Avr II</i>	BsrBR I	<i>Eco81 I</i>	<i>Pac I</i>	SnaB I	
Bal I	<i>BsrG I</i>	<i>EcoN I</i>	<i>PflM I</i>	<i>SpI I</i>	
BamH I	BssH II	EcoR I	<i>PinA I</i>	<i>Srf I</i>	
<i>Bbe I</i>	<i>Bst1107 I</i>	<i>Ehe I</i>	<i>Pme I</i>	Stu I	
<i>BbrP I</i>	Bst98 I	<i>Fse I</i>	<i>Pml I</i>	<i>Swa I</i>	
<i>Bbs I</i>	BstE II	Hind III	<i>PpuM I</i>	Tth111 I	

Table 3. Restriction Enzymes That Cut the pGEM®-5Zf(–) Vector 6 or More Times.

<i>Aci I</i>	<i>BstU I</i>	Hinf I	<i>Mnl I</i>	Sau3A I
Alu I	Cfo I	Hpa II	<i>Mse I</i>	Sau96 I
<i>Bbv I</i>	Dpn I	<i>Hph I</i>	Msp I	<i>ScrF I</i>
BsaO I	<i>Dpn II</i>	Hsp92 II	MspA1 I	<i>SfaN I</i>
Bsp1286 I	<i>Eae I</i>	<i>Mae II</i>	Nde II	Tru9 I
<i>Bsr I</i>	<i>Fnu4H I</i>	<i>Mae III</i>	<i>Nla III</i>	Xho II
BsrS I	Hae III	Mbo I	<i>Nla IV</i>	
Bst71 I	Hha I	Mbo II	<i>Ple I</i>	

Note: The enzymes listed in boldface type are available from Promega.

V. pGEM®-5Zf(–) Vector Sequence

The sequence shown corresponds to RNA synthesized by T7 RNA polymerase and is complementary to RNA synthesized by SP6 RNA polymerase. The strand shown is the same as the ssDNA strand produced by this vector.

```

1  GGGCGAATTG  GGCCCGACGT  CGCATGCTCC  CGGCCGCCAT  GGCCGCGGGA
51  TATCACTAGT  GCGGCCGCCT  GCAGGTCGAC  CATATGGGAG  AGCTCCCAAC
101 GCGTTGGATG  CATAGCTTGA  GTATTCTATA  GTGTCACCTA  AATAGCTTGG
151 CGTAATCATG  GTCATAGCTG  TTTCTGTGT  GAAATTGTTA  TCCGCTCACA
201 ATTCCACACA  ACATACGAGC  CGGAAGCATA  AAGTGTAAG  CCTGGGGTGC
251 CTAATGAGTG  AGCTAACTCA  CATTAATTGC  GTTGCGCTCA  CTGCCCCTT
301 TCCAGTCGGG  AAACCTGTCG  TGCCAGCTGC  ATTAATGAAT  CGGCCAACGC
351 GCGGGGAGAG  GCGGTTTGCG  TATTGGGCGC  TCTTCCGCTT  CCTCGCTCAC
401 TGA CTGCTG  CGCTCGGTCTG  TTCGGCTGCG  GCGAGCGGTA  TCAGCTCACT
451 CAAAGGCGGT  AATACGGTTA  TCCACAGAAT  CAGGGGATAA  CGCAGGAAAG
501 AACATGTGAG  CAAAAGGCCA  GCAAAAGGCC  AGGAACCGTA  AAAAGGCCGC
551 GTTGCTGGCG  TTTTTCATA  GGCTCCGCC  CCCTGACGAG  CATCACAAAA
601 ATCGACGCTC  AAGTCAGAGG  TGGCGAAACC  CGACAGGACT  ATAAAGATAC

```

pGEM®-5Zf(-) Vector Sequence (continued)

651	CAGGCGTTTC	CCCCTGGAAG	CTCCCTCGTG	CGCTCTCCTG	TTCCGACCCT
701	GCCGCTTACC	GGATACCTGT	CCGCCTTTCT	CCCTTCGGGA	AGCGTGGCGC
751	TTTCTCATAG	CTCACGCTGT	AGGTATCTCA	GTTCCGGTGTA	GGTCGTTTCG
801	TCCAAGCTGG	GCTGTGTGCA	CGAACCCCCC	GTTCAGCCCCG	ACCGCTGCGC
851	CTTATCCGGT	AACTATCGTC	TTGAGTCCAA	CCCGGTAAGA	CACGACTTAT
901	CGCCACTGGC	AGCAGCCACT	GGTAACAGGA	TTAGCAGAGC	GAGGTATGTA
951	GGCGGTGCTA	CAGAGTTCTT	GAAGTGGTGG	CCTAACTACG	GCTACACTAG
1001	AAGAACAGTA	TTTGGTATCT	GCGCTCTGCT	GAAGCCAGTT	ACCTTCGGAA
1051	AAAGAGTTGG	TAGCTCTTGA	TCCGGCAAAC	AAACCACCGC	TGGTAGCGGT
1101	GGTTTTTTTTG	TTTGCAAGCA	GCAGATTACG	CGCAGAAAAA	AAGGATCTCA
1151	AGAAGATCCT	TTGATCTTTT	CTACGGGGTC	TGACGCTCAG	TGGAACGAAA
1201	ACTCACGTTA	AGGGATTTTG	GTCATGAGAT	TATCAAAAAG	GATCTTCACC
1251	TAGATCCTTT	TAAATTAAAA	ATGAAGTTTT	AAATCAATCT	AAAGTATATA
1301	TGAGTAAACT	TGGTCTGACA	GTTACCAATG	CTTAATCAGT	GAGGCACCTA
1351	TCTCAGCGAT	CTGTCTATTT	CGTTCATCCA	TAGTTGCCTG	ACTCCCCGTC
1401	GTGTAGATAA	CTACGATACG	GGAGGGCTTA	CCATCTGGCC	CCAGTGCTGC
1451	AATGATACCG	CGAGACCCAC	GCTCACCGGC	TCCAGATTTA	TCAGCAATAA
1501	ACCAGCCAGC	CGGAAGGGCC	GAGCGCAGAA	GTGGTCCTGC	AACTTTATCC
1551	GCCTCCATCC	AGTCTATTAA	TTGTTGCCGG	GAAGCTAGAG	TAAGTAGTTC
1601	GCCAGTTAAT	AGTTTGCGCA	ACGTTGTTGC	CATTGCTACA	GGCATCGTGG
1651	TGTCACGCTC	GTCGTTTGGT	ATGGCTTCAT	TCAGCTCCGG	TTCCCAACGA
1701	TCAAGGCGAG	TTACATGATC	CCCCATGTTG	TGCAAAAAG	CGGTTAGCTC
1751	CTTCGGTCCT	CCGATCGTTG	TCAGAAGTAA	GTTGGCCGCA	GTGTTATCAC
1801	TCATGGTTAT	GGCAGCACTG	CATAATTCTC	TTACTGTCAT	GCCATCCGTA
1851	AGATGCTTTT	CTGTGACTGG	TGAGTACTCA	ACCAAGTCAT	TCTGAGAATA
1901	GTGTATGCGG	CGACCGAGTT	GCTCTTGCCC	GGCGTCAATA	CGGGATAATA
1951	CCGCGCCACA	TAGCAGAACT	TTAAAAGTGC	TCATCATTTG	AAAACGTTCT
2001	TCGGGGCGAA	AACTCTCAAG	GATCTTACCG	CTGTTGAGAT	CCAGTTCGAT
2051	GTAACCCACT	CGTGACCCCA	ACTGATCTTC	AGCATCTTTT	ACTTTCACCA
2101	GCGTTTCTGG	GTGAGCAAAA	ACAGGAAGGC	AAAATGCCGC	AAAAAAGGGA
2151	ATAAGGGCGA	CACGGAAATG	TTGAATACTC	ATACTCTTCC	TTTTTCAATA
2201	TTATTGAAGC	ATTTATCAGG	GTTATTGTCT	CATGAGCGGA	TACATATTTG
2251	AATGTATTTA	GAAAAATAAA	CAAATAGGGG	TTCCGCGCAC	ATTTCCCCGA
2301	AAAGTGCCAC	CTGTATGCGG	TGTGAAATAC	CGCACAGATG	CGTAAGGAGA
2351	AAATACCGCA	TCAGGACGCG	CCCTGTAGCG	GCGCATTAAG	CGCGGCGGGT
2401	GTGGTGTTTA	CGCGCAGCGT	GACCGCTACA	CTTGCCAGCG	CCCTAGCGCC
2451	CGCTCCTTTC	GCTTTCTTCC	CTTCCTTTCT	CGCCACGTTC	GCCGGCTTTC
2501	CCCGTCAAGC	TCTAAATCGG	GGGCTCCCTT	TAGGGTTCCG	ATTTAGTGCT

pGEM®-5Zf(-) Vector Sequence (continued)

```

2551  TTACGGCACC  TCGACCCCAA  AAAACTTGAT  TAGGGTGATG  GTTCACGTAG
2601  TGGGCCATCG  CCCTGATAGA  CGGTTTTTCG  CCCTTTGACG  TTGGAGTCCA
2651  CGTTCTTTAA  TAGTGGACTC  TTGTTCCAAA  CTGGAACAAC  ACTCAACCCT
2701  ATCTCGGTCT  ATTCTTTTGA  TTTATAAGGG  ATTTTGCCGA  TTTCGGCCTA
2751  TTGGTTAAAA  AATGAGCTGA  TTTAACAAAA  ATTTAACGCG  AATTTTAACA
2801  AAATATTAAC  GCTTACAATT  TCCATTCGCC  ATTCAGGCTG  CGCAACTGTT
2851  GGGAAGGGCG  ATCGGTGCGG  GCCTCTTCGC  TATTACGCCA  GCTGGCGAAA
2901  GGGGGATGTG  CTGCAAGGCG  ATTAAGTTGG  GTAACGCCAG  GGTTTTCCCA
2951  GTCACGACGT  TGTAAAACGA  CGGCCAGTGA  ATTGTAATAC  GACTCACTAT
3001  A

```

VI. Related Products

pGEM® Vectors

Product	Size	Cat.#
pGEM®-3Z Vector(a)	20µg	P2151
pGEM®-4Z Vector(a)	20µg	P2161
pGEM®-3Zf(+) Vector(a)	20µg	P2271
pGEM®-3Zf(-) Vector(a)	20µg	P2261
pGEM®-5Zf(+) Vector(a)	20µg	P2241
pGEM®-7Zf(+) Vector(a)	20µg	P2251
pGEM®-7Zf(-) Vector(a)	20µg	P2371
pGEM®-9Zf(-) Vector(a)	20µg	P2391
pGEM®-11Zf(+) Vector(a)	20µg	P2411
pGEM®-11Zf(-) Vector(a)	20µg	P2421
pGEM®-13Zf(+) Vector(a)	20µg	P2541

All pGEM®-Zf Vectors are provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain the vector and are not competent.

Other Vectors

Product	Size	Cat.#
pSP64 Poly(A) Vector	20µg	P1241
pSP72 Vector(a)	20µg	P2191
pSP73 Vector(a)	20µg	P2221

Sequencing Primers

Product	Size	Cat.#
SP6 Promoter Primer	2µg	Q5011
T7 Promoter Primer	2µg	Q5021
pUC/M13 Primer, Reverse (17mer)	2µg	Q5401
pUC/M13 Primer, Forward (17mer)	2µg	Q5391
pUC/M13 Primer, Forward (24mer)	2µg	Q5601
pUC/M13 Primer, Reverse (22mer)	2µg	Q5421



Related Products (continued)

Erase-a-Base® Systems

Product	Cat.#
Erase-a-Base® System (plus vectors)	E5850
Erase-a-Base® System (minus vectors and bacterial strain)	E5750

Riboprobe® in vitro Transcription Systems

Product	Cat.#
Riboprobe® System - SP6(b,c)	P1420
Riboprobe® System - T3(b,c)	P1430
Riboprobe® System - T7(b,c)	P1440

TnT® Quick Coupled Transcription/Translation Systems

Product	Cat.#
TnT® T7 Quick Coupled Transcription/Translation System(b,c,d,e,f)	L1170
TnT® T7 Quick Coupled Transcription/Translation System, Trial Size(b,c,d,e,f)	L1171
TnT® SP6 Quick Coupled Transcription/Translation System(b,c,d,e,f)	L2080
TnT® SP6 Quick Coupled Transcription/Translation System, Trial Size(b,c,d,e,f)	L2081

VII. Reference

1. Yanish-Perron, C. *et al.* (1985) Improved M13 phage cloning vectors and host strains: nucleotide sequences of the M13mp18 and pUC19 vectors. *Gene* **33**, 103.

(a)U.S. Pat. No. 4,766,072 has been issued to Promega Corporation for transcription vectors having two different bacteriophage RNA polymerase promoter sequences separated by a series of unique restriction sites into which foreign DNA can be inserted.

(b)U.S. Pat. No. 5,552,302, European Pat. No. 0 422 217 and Australian Pat. No. 646803 have been issued to Promega Corporation for the methods and compositions for production of human recombinant placental ribonuclease inhibitor.

(c)U.S. Pat. Nos. 4,966,964, 5,019,556 and 5,266,687, which claim vectors encoding a portion of human placental ribonuclease inhibitor, are exclusively licensed to Promega Corporation.

(d)U.S. Pat. Nos. 5,283,179, 5,641,641, 5,650,289, Australian Pat. No. 649289 and European Pat. No. 0 553 234 have been issued to Promega Corporation for a firefly luciferase assay method, which affords greater light output with improved kinetics as compared to the conventional assay. Certain applications of this product may require licenses from others.

(e)U.S. Pat. Nos. 5,324,637, 5,492,817 and 5,665,563, European Pat. No. 0 566 714 B1, Australian Pat. No. 660329 and Japanese Pat. No. 2904583 have been issued to Promega Corporation for coupled transcription/translation systems that use RNA polymerases and eukaryotic lysates.

(f)The method of recombinant expression of *Coleoptera* luciferase is covered by U.S. Pat. Nos. 5,583,024, 5,674,713 and 5,700,673.

(g)U.S. Pat. No. 5,391,487 has been issued to Promega Corporation for Restriction Endonuclease *Sgf I*.

©1988–2000 Promega Corporation. All Rights Reserved.

Erase-a-Base, pGEM, Riboprobe and TnT are trademarks of Promega Corporation and are registered with the U.S. Patent and Trademark Office.

DNASTAR is a registered trademark of DNASTAR, Inc. GenBank is a registered trademark of U.S. Department of Health and Human Services.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.



Promega Corporation

2800 Woods Hollow Road	
Madison, WI 53711-5399	USA
Telephone	608-274-4330
Fax	608-277-2516
Internet	www.promega.com